






APRIL 2018


II INFORMATION FOR LABORATORIES AND PRACTITIONERS

M1.4 SCREW GUIDE FOR CUSTOMISED SIMEDA® AND CONNECT+® MULTIPLE PROSTHESES WITH ANGULATED ACCESS

- Screw for screw-retained prostheses on Multi-Unit Axiom® BL - **Simeda®** abutments

Permanent screw				Laboratory screw	
Zirconia prosthesis Sina Z®/T®/ML NEW		Titanium or CoCr prosthesis		Titanium, CoCr or zirconia Sina Z®/T®/ML prosthesis	
Platform 4.0	Platform 4.8	Platform 4.0	Platform 4.8	Platform 4.0	Platform 4.8
M1.4 titanium screw  Ref. MUAA141		M1.4 bicolour screw  Ref. MUAA140		M1.4 blue screw (x 4)  Ref. MUAA142-4	

- Screw for screw-retained prostheses on Multi-Unit NobelBiocare® - **Connect+®** abutments

Permanent screw		Laboratory screw	
For Titanium, CoCr or zirconia Sina Z®/T®/ML prosthesis NEW		Titanium, CoCr or zirconia Sina Z®/T®/ML prosthesis	
RP Platform	WP Platform	RP Platform	WP Platform
M1.4 titanium screw  Ref. MUAA141	M1.8 titanium screw  Ref. MUAA181	M1.4 blue screw (x 4)  Ref. MUAA142-4	M1.8 blue screw (x 1)  Ref. MUAA182